

Human GOLPH2/GOLM1 protein (Recombinant) (C-His) (STJP001499)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51280
Gene Symbol	GOLM1
UniProt ID	GOLM1_HUMAN
Immunogen Region	Ser36-Leu401
Immunogen Sequence	SSRSVDLQTRIMELEGRVRR AAAERGAVELKKNEFQGELE KQREQLDKIQSSSHNFQLESV NKLYQDEKAVLVNNTTGER LIRVLQDQLKTLQRNYGRLQ QDVLQFQKNQTNLERKFSYD LSQCINQMKEVKEQCCEERIE EVTCKGNEAVASRDLSENNDRQQLQALSEPQPRQLQAAGL PHTEVPQGKGNVLGNSKSQT PAPSSEVVLDSKRQVEKEET NEIQVVNEEPQRDRLPQEP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

