

Mouse CD69 protein (Recombinant) (C-His) (STJP001503)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	12515
Gene Symbol	Cd69
UniProt ID	CD69_MOUSE
Immunogen Region	Asn62-Arg199
Immunogen Sequence	NVGKYNCPGLYEKLESSDHH VATCKNEWISYKRTCYFFST TTKSWALAQRSCSEDAATLA VIDSEKDMTFLKRYSGELEH WIGLKNEANQTWKWANGKEF NSWFNLTGSGRCVSVNHKNV TAVDCEANFHWVCSKPSR

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

