

## SARS-CoV-2 Nucleocapsid (B.1.617.2/Delta) protein (Recombinant) (C-His&Avi) (STJP001504)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20 mM Tris, 150 mM NaCl, 10% Glycerol, pH 8.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Immunogen Region	Ser2-Ala419
Immunogen Sequence	SDNGPQNQRNAPRITFGGPS DSTGSNQNGERSGARSKQRR PQGLPNNTASWFTALTQHGK EGLKFPRGQGVPIINTNSSPD DQIGYYRRATRIRGGDGKM KDLSPRWYFYLLGTGPEAGL PYGANKDGIWVATEGALNT PKDHIGTRNPANNAIVLQL PQGTTLPKGFYAEGRGGGSSQ ASSRSSRSRNSSRNSTPGS SMGTSPARMAGNGGDAALAL LLLDRLNQLESKMSGKGQQ

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

