

Human IFNAR1 protein (Recombinant) (C-His) (STJP001505)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3454
Gene Symbol	IFNAR1
UniProt ID	INAR1_HUMAN
Immunogen Region	Lys28-Lys436
Immunogen Sequence	KNLKSPQKVEVDIIDDNFIL RWNRSDESVGNVTFSDYQK TGMDNWIKLSGCQNITSTKC NFSSLKLNVEEIKLRIRAE KENTSSWYEVDSFTPFRAQ IGPPEVHLEAEDKAIVIHIS PGTKDSVMWALDGLSFTYSL VIWKNSSGVEERIEIYSRH KIYKLPETTYCLKVKAALL TSWKIGVYSPVHCIKTTVEN ELPPPENIEVSVQNQNYVLK WDYTYANMTFQVQWLHAFLL

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

