

Mouse CD8a protein (Recombinant) (C-hFc) (STJP001508)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	12525
Gene Symbol	Cd8a
UniProt ID	CD8A_MOUSE
Immunogen Region	Lys28-Tyr196
Immunogen Sequence	KPQAPELRIFPKMDAELGQ KVDLVCEVLGVSQGC SWLF QNSSSKLPQPTFVVYMASSH NKITWDEKLNSSKLF SAMRD TNNKYVLT LNKFSKENEGYY FCSVISNSVMYFSSVVPVLQ KVNSTTTKPVL RTPSPVHPT GTSQPQRPEDCRPRGSVKGT GLDFACDIY

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

