

Human KLRG1/CLEC15A/CD369 protein (Recombinant) (N-raFc) (STJP001515)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10219
Gene Symbol	KLRG1
UniProt ID	KLRG1_HUMAN
Immunogen Region	Leu60-Phe195
Immunogen Sequence	LCQGSNYSTCASCPCPDWR MKYGNHCYFVSVEEKDWNSS LEFCLARDSHLLVITDNQEM SLLQVFLSEAFWIGLRNNS GWRWEDGSPLNFSRISSNSF VQTCGAINKNGLQASSCEVP LHWVCKKCPFADQALF

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

