

Human Syndecan-1/SDC1/CD138 protein (Recombinant) (C-raFc) (STJP001524)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6382
Gene Symbol	SDC1
UniProt ID	SDC1_HUMAN
Immunogen Region	Gln23-Glu251
Immunogen Sequence	QIVATNLPPEDQDGSDDSD NFGSGAGALQDITLSQQTP STWKDTQLLTAIPTSPEPTG LEATAASTSTLPAGEGPKEG EAVVLPEVEPGLTAREQEAT PRPRETTQLPTTHLASTTTA TTAQEPATSHPHRDMQPGHH ETSTPAGPSQADLHTPHTED GGPSATERAAEDGASSQLPA AEGSGEQDFTFETSGENTAV VAVEPDRRNQSPVDQGATGA SQGLLDRKE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

