

Human Lipocalin-1/LCN1 protein (Recombinant) (C-His) (STJP001530)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3933
Gene Symbol	LCN1
UniProt ID	LCN1_HUMAN
Immunogen Region	His19-Asp176
Immunogen Sequence	HHLLASDEEIQDVSGTWYLK AMTVDREFPEMNLESVTPMT LTTLEGGNLEAKVTMLISGR CQEVKAVLEKTDEPGKYTAD GGKHVAYIIRSHVKDHYIFY CEGELHGKPVRGVKLVGRDP KNNLEALEDFEKAAGARGLS TESILIPRQSETCSPGSD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

