

Human MCP-1/CCL2 protein (Recombinant) (C-His-Avi) {Biotin} (STJP001547)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES



Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4.
Molecular Weight	Recombinant Human MCP-1/CCL2 is a protein consisting of amino acid residues [Gln24-Thr99]. The protein has a predicted MW of 11.5 kDa. Due to glycosylation, the protein migrates to 15-20 kDa based on Tris-Bis PAGE result..
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.

TARGET INFORMATION

Gene ID	6347
Gene Symbol	CCL2
UniProt ID	CCL2_HUMAN
Immunogen Region	Gln24-Thr99

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Biotinylated Human MCP-1 on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.

The purity of Biotinylated Human MCP-1 is greater than 95% as determined by SEC-HPLC.