

Human CD30/TNFRSF8 protein (Recombinant) (C-His-Avi) {Biotin} (STJP001550)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4.
Molecular Weight	Recombinant Human CD30/TNFRSF8 is a protein consisting of amino acid residues [Phe19-Lys379]. The protein has a predicted MW of 41.3 kDa. Due to glycosylation, the protein migrates to 70-100 kDa based on Tris-Bis PAGE result..
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.

TARGET INFORMATION


Gene ID	943
Gene Symbol	TNFRSF8
UniProt ID	TNR8_HUMAN
Immunogen Region	Phe19-Lys379

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Biotinylated CD30 on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.



The purity of Biotinylated CD30 is greater than 95% as determined by SEC-HPLC.



Immobilized Anti-CD30 Antibody, hFc Tag at 0.5 μ g/ml (100 μ g/ml) on the plate. Dose response curve for Biotinylated Human CD30, His Tag with the EC50 of 49 ng/ml determined by ELISA.