

Human Notch 3 protein (Recombinant) (C-His-Avi) {Biotin} (STJP001557)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES



Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4.
Molecular Weight	Recombinant Human Notch 3 is a protein consisting of amino acid residues [Ala40-Glu467]. The protein has a predicted MW of 47.7 kDa. Due to glycosylation, the protein migrates to 55-65 kDa based on Tris-Bis PAGE result..
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.

TARGET INFORMATION

Gene ID	4854
Gene Symbol	NOTCH3
UniProt ID	NOTC3_HUMAN
Immunogen Region	Ala40-Glu467

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Biotinylated Human Notch 3 on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.

The purity of Biotinylated Human Notch 3 is greater than 95% as determined by SEC-HPLC.