

Mouse ALCAM/CD166 protein (Recombinant) (C-His-Avi) {Biotin} (STJP001569)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4.
Molecular Weight	Recombinant Mouse ALCAM/CD166 is a protein consisting of amino acid residues [Tyr28-Lys527]. The protein has a predicted MW of 58.9 kDa. Due to glycosylation, the protein migrates to 75-85 kDa based on Tris-Bis PAGE result..
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.

TARGET INFORMATION


Gene ID	11658
Gene Symbol	Alcam
UniProt ID	CD166_MOUSE
Immunogen Region	Tyr28-Lys527

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--



Biotinylated Mouse ALCAM on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.



The purity of Biotinylated Mouse ALCAM is greater than 95% as determined by SEC-HPLC.