

Human FcRH5/FcRL5 protein (Recombinant) (C-His-Avi) {Biotin} (STJP001583)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4.
Molecular Weight	Recombinant Human FcRH5/FcRL5 is a protein consisting of amino acid residues [Gln16-Gly851]. The protein has a predicted MW of 94.2 kDa. Due to glycosylation, the protein migrates to 110-120 kDa based on Tris-Bis PAGE result..
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.

TARGET INFORMATION


Gene ID	83416
Gene Symbol	FCRL5
UniProt ID	FCRL5_HUMAN
Immunogen Region	Gln16-Gly851

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Biotinylated Human FcRH5 on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.



Immobilized Anti-FcRH5 Antibody at 2 μ g/ml (100 μ g/ml) on the plate. Dose response curve for Biotinylated Human FcRH5, His Tag with the EC50 of 0.23 μ g/ml determined by ELISA.