

Human SIRP Beta 1 Isoform 3 protein (Recombinant) (C-His-Avi) {Biotin} (STJP001599)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4.
Molecular Weight	Recombinant Human SIRP Beta 1 Isoform 3 is a protein consisting of amino acid residues [Glu30-Leu371], The protein has a predicted MW of 39.9 kDa. Due to glycosylation, the protein migrates to 55-68 kDa based on Tris-Bis PAGE result..
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.

TARGET INFORMATION


Gene ID	10326
Gene Symbol	SIRPB1
UniProt ID	SIRBL_HUMAN
Immunogen Region	Glu30-Leu371

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Biotinylated Human SIRP Beta 1 Isoform 3 on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.



The purity of Biotinylated Human SIRP Beta 1 Isoform 3 is greater than 95% as determined by SEC-HPLC.