

Human LILRA4/CD85g protein (Recombinant) (C-His-Avi) {Biotin} (STJP001600)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4.
Molecular Weight	Recombinant Human LILRA4/CD85g is a protein consisting of amino acid residues [Glu24-Asn446]. The protein has a predicted MW of 49.5 kDa. Due to glycosylation, the protein migrates to 60-70 kDa based on Tris-Bis PAGE result..
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.

TARGET INFORMATION


Gene ID	23547
Gene Symbol	LILRA4
UniProt ID	LIRA4_HUMAN
Immunogen Region	Glu24-Asn446

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--



Biotinylated Human LILRA4 on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.



The purity of Biotinylated Human LILRA4 is greater than 95% as determined by SEC-HPLC.