

Human ENPP-3 protein (Recombinant) (N-His-Avi) {Biotin} (STJP001611)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES



Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4.
Molecular Weight	Recombinant Human ENPP-3 is a protein consisting of amino acid residues [Leu48-Ile875]. The protein has a predicted MW of 99.97 kDa. Due to glycosylation, the protein migrates to 110-130 kDa based on Tris-Bis PAGE result..
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.

TARGET INFORMATION

Gene ID	5169
Gene Symbol	ENPP3
UniProt ID	ENPP3_HUMAN
Immunogen Region	Leu48-Ile875

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Biotinylated Human ENPP-3 on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.

The purity of Biotinylated Human ENPP-3 is greater than 95% as determined by SEC-HPLC.