

Human KIR2DL5 protein (Recombinant) (C-His-Avi) {Biotin} (STJP001615)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4.
Molecular Weight	Recombinant Human KIR2DL5 is a protein consisting of amino acid residues [His22-His240], The protein has a predicted MW of 26.3 kDa. Due to glycosylation, the protein migrates to 45-52 kDa based on Tris-Bis PAGE result..
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.

TARGET INFORMATION

Immunogen Region	His22-His240
------------------	--------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--



Biotinylated Human KIR2DL5 on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.