

Human TGFBR1 protein (Recombinant) (C-mFc-Avi) {Biotin} (STJP001628)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4.
Molecular Weight	Recombinant Human TGFBR1 is a protein consisting of amino acid residues [Leu34-Glu125]. The protein has a predicted MW of 38.2 kDa. Due to glycosylation, the protein migrates to 50-55 kDa based on Tris-Bis PAGE result..
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.

TARGET INFORMATION

Gene ID	7046
Gene Symbol	TGFBRI1
UniProt ID	TGFR1_HUMAN
Immunogen Region	Leu34-Glu125

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--



Biotinylated Human TGFBR1 on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.