

## Human TFPI protein (Recombinant) (C-His-Avi) {Biotin} (STJP001639)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

### PRODUCT PROPERTIES


Formulation	Lyophilised from 0.2 $\mu$ m filtered PBS solution pH7.4.
Molecular Weight	Recombinant Human TFPI is a protein consisting of amino acid residues [Asp29-Lys282]. The protein has a predicted MW of 32.1 kDa. Due to glycosylation, the protein migrates to 40-50 kDa based on Tris-Bis PAGE result..
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.

### TARGET INFORMATION

Gene ID	<a href="#">7035</a>
Gene Symbol	<a href="#">TFPI</a>
UniProt ID	<a href="#">TFPI1_HUMAN</a>
Immunogen Region	Asp29-Lys282

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--



Biotinylated Human TFPI on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.