

Human BAFFR/TNFRSF13C protein (Recombinant) (C-His-Avi) {Biotin} (STJP001656)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4.
Molecular Weight	Recombinant Human BAFFR/TNFRSF13C is a protein consisting of amino acid residues [Ser7-Ala71]. The protein has a predicted MW of 9.5 kDa. Due to glycosylation, the protein migrates to 20-25 kDa based on Tris-Bis PAGE result..
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.

TARGET INFORMATION


Gene ID	115650
Gene Symbol	TNFRSF13C
UniProt ID	TR13C_HUMAN
Immunogen Region	Ser7-Ala71

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Biotinylated Human BAFFR on Tris-Bis PAGE under reduced conditions. The purity is greater than 95%.



The purity of Biotinylated Human BAFFR is greater than 95% as determined by SEC-HPLC.



Immobilized Anti-BAFFR Antibody, hFc Tag at 0.5 μ g/ml (100 μ g/ml) on the plate. Dose response curve for Biotinylated Human BAFFR, His Tag with the EC50 of 3.4 μ g/ml determined by ELISA.