

Human LILRA2/CD85h/ILT1 protein (Recombinant) (C-His-Avi) {Biotin} (STJP001667)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES



Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4.
Molecular Weight	Recombinant Human LILRA2/CD85h/ILT1 is a protein consisting of amino acid residues [Gly24-Asn449], The protein has a predicted MW of 49.9 kDa. Due to glycosylation, the protein migrates to 70-80 kDa based on Tris-Bis PAGE result..
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.

TARGET INFORMATION

Gene ID	11027
Gene Symbol	LILRA2
UniProt ID	LIRA2_HUMAN
Immunogen Region	Gly24-Asn449

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--



Biotinylated Human LILRA2 on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.

The purity of Biotinylated Human LILRA2 is greater than 95% as determined by SEC-HPLC.