

Human LILRA1/CD85i/LIR-6 protein (Recombinant) (C-His-Avi) {Biotin} (STJP001689)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4.
Molecular Weight	Recombinant Human LILRA1/CD85i/LIR-6 is a protein consisting of amino acid residues [Pro17-Asn461]. The protein has a predicted MW of 51.4 kDa. Due to glycosylation, the protein migrates to 70-80 kDa based on Tris-Bis PAGE result..
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.

TARGET INFORMATION


Gene ID	11024
Gene Symbol	LILRA1
UniProt ID	LIRA1_HUMAN
Immunogen Region	Pro17-Asn461

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Biotinylated Human LILRA1 on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.



The purity of Biotinylated Human LILRA1 is greater than 95% as determined by SEC-HPLC.