

Human BCMA/TNFRSF17 protein (Recombinant) (C-His-Avi) {Biotin} (STJP001691)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4.
Molecular Weight	Recombinant Human BCMA/TNFRSF17 is a protein consisting of amino acid residues [Met1-Ala54]. The protein has a predicted MW of 8.9 kDa. Due to glycosylation, the protein migrates to 12 kDa and 15-17 kDa based on Tris-Bis PAGE result..
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.

TARGET INFORMATION


Gene ID	608
Gene Symbol	TNFRSF17
UniProt ID	TNR17_HUMAN
Immunogen Region	Met1-Ala54

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Biotinylated Human BCMA on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.



The purity of Biotinylated Human BCMA is greater than 95% as determined by SEC-HPLC.



Immobilized Anti-BCMA Antibody, hFc Tag at 0.2 μ g/ml (100 μ g I/Well) on the plate. Dose response curve for Biotinylated Human BCMA, His Tag with the EC50 of 1.8ng/ml determined by ELISA.