

Cynomolgus BAFF/TNFSF13B/CD257 Trimer protein (Recombinant) (N-His-Avi, N-Flag) {Biotin} (STJP001729)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4.
Molecular Weight	Recombinant Cynomolgus BAFF/TNFSF13B/CD257 Trimer is a protein consisting of amino acid residues [Thr141-Leu285]. The protein has a predicted MW of 54.2 kDa. Due to glycosylation, the protein migrates to 55-60 kDa based on Tris-Bis PAGE result..
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.

TARGET INFORMATION


Immunogen Region	Thr141-Leu285
------------------	---------------

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Biotinylated Cynomolgus BAFF Trimer on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.



The purity of Biotinylated Cynomolgus BAFF Trimer is greater than 95% as determined by SEC-HPLC.



Immobilized Biotinylated Cynomolgus BAFF Trimer, His Tag at 1 Mu g/ml (100 Mu I/Well) on the plate. Dose response curve for Human BAFFR, hFc Tag with the EC50 of 0.39 Mu g/ml determined by ELISA.