

Human ANGPTL3/Angiopoietin-like 3 protein (Recombinant) (C-His-Avi) {Biotin} (STJP001735)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4.
Molecular Weight	Recombinant Human ANGPTL3/Angiopoietin-like 3 is a protein consisting of amino acid residues [Ser17-Lys219]. The protein has a predicted MW of 26.6 kDa. Due to glycosylation, the protein migrates to 30-38 kDa based on Tris-Bis PAGE result..
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.

TARGET INFORMATION


Gene ID	27329
Gene Symbol	ANGPTL3
UniProt ID	ANGL3_HUMAN
Immunogen Region	Ser17-Lys219

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Biotinylated Human ANGPTL3 on Tris-Bis PAGE under reduced conditions. The purity is greater than 95%.



Immobilized Anti-ANGPTL3 Antibody, hFc Tag at 1 μ g/ml (100 μ g/ml) on the plate. Dose response curve for Biotinylated Human ANGPTL3, His Tag with the EC50 of 5.0 ng/ml determined by ELISA.