

Mouse MSLN/Mesothelin protein (Recombinant) (C-His) {Biotin} (STJP001749)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES



Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4.
Molecular Weight	Recombinant Mouse MSLN/Mesothelin (Primary Amine Labeling) is a protein consisting of amino acid residues [Asp298-Ser600]. The protein has a predicted MW of 35.2 kDa. Due to glycosylation, the protein migrates to 45-60 kDa based on Tris-Bis PAGE resu
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.

TARGET INFORMATION

Gene ID	56047
Gene Symbol	Msln
UniProt ID	MSLN_MOUSE
Immunogen Region	Asp298-Ser600

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Biotinylated Mouse MSLN on Tris-Bis PAGE under reduced condition. The purity is greater than 95% as determined by SEC-HPLC.