

Human CD155/PVR protein (Recombinant) (C-His-Avi) {Biotin} (STJP001755)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4.
Molecular Weight	Recombinant Human CD155/PVR is a protein consisting of amino acid residues [Trp21-Asn343]. The protein has a predicted MW of 38 kDa. Due to glycosylation, the protein migrates to 55-70 kDa based on Tris-Bis PAGE result..
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.

TARGET INFORMATION


Gene ID	5817
Gene Symbol	PVR
UniProt ID	PVR_HUMAN
Immunogen Region	Trp21-Asn343

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Biotinylated Human CD155 on Tris-Bis PAGE under reduced conditions. The purity is greater than 95%.



The purity of Biotinylated Human CD155 is greater than 95% as determined by SEC-HPLC.



Immobilized Biotinylated Human CD155, His Tag at 1 μ g/ml (100 μ g/ml) on the streptavidin precoated plate (5 μ g/ml). Dose response curve for Human TIGIT, hFc Tag with the EC₅₀ of 4.3ng/ml determined by ELISA.