

HLA-A*02:01 gp100 peptide complex Tetramer peptide (Recombinant) (STJP001790)

GENERAL INFORMATION

Product Type	Biomolecules
Host / Source	HEK293 cells NA

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4 NA
Molecular Weight	Recombinant HLA-A/GP100 Tetramer has a calculated mass of approximately 50kDa. Due to glycosylation migrates as an approximately 60 kDa protein under reducing conditions in SDS-PAGE. NA
Storage Instruction	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C C. NA

TARGET INFORMATION

Immunogen Region	Human HLA-A*02:01 AA25-305/B2M protein AA2-119
------------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--