

Human 2B4/CD244/SLAMF4 Protein protein (Recombinant) (STJP001802)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells NA

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 μ m filtered PBS solution pH7.4
Molecular Weight	Recombinant Human 2B4/CD244/SLAMF4 is a protein consisting of amino acid residues [Cys22-Arg221]. The protein has a predicted MW of 49 kDa. Due to glycosylation, the protein migrates to 68-75 kDa based on Tris-Bis PAGE result..
Storage Instruction	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C C.

TARGET INFORMATION

Gene ID	51744
Gene Symbol	CD244
UniProt ID	CD244_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--