

## Human Cortical thymocyte-like protein protein (Recombinant) (STJP001855)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

### PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Molecular Weight	Recombinant Human Cortical thymocyte-like protein (VSIG2) is a protein consisting of 227 amino acid residue subunits, due to glycosylation migrates as an approximately as 35kDa band protein on reduced SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at-20°C to-70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

### TARGET INFORMATION

Gene ID	<a href="#">23584</a>
Gene Symbol	<a href="#">VSIG2</a>
UniProt ID	<a href="#">VSIG2_HUMAN</a>
Immunogen Region	VAL24-ALA243

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--