

## Human uPAR protein (Recombinant) (STJP001869)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

### PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Molecular Weight	Recombinant human Urokinase-type plasminogen activator receptor (uPAR) is a protein consisting of 390 amino acid residue subunits, due to glycosylation migrates as an approximately 50kDa bands protein on reduced SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

### TARGET INFORMATION

Gene ID	<a href="#">5329</a>
Gene Symbol	<a href="#">PLAUR</a>
UniProt ID	<a href="#">UPAR_HUMAN</a>
Immunogen Region	LEU23-ARG303

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--