

Human Urokinase-type plasminogen activator protein (Recombinant) (STJP001870)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 M filtered PBS solution, pH7.2, 5% Trehalose.
Molecular Weight	Recombinant human Urokinase-type plasminogen activator (uPA) is a protein consisting of 420 amino acid residue subunits, due to glycosylation and cleavage migrates as an approximately 18kDa, 32kDa, 50kDa bands protein on reduced SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	5328
Gene Symbol	PLAU
UniProt ID	UROK_HUMAN
Immunogen Region	SER21-LEU431

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--