

Mouse c-MET protein (Recombinant) (STJP001877)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 M filtered PBS solution, pH7.2, 5% Trehalose.
Molecular Weight	Recombinant mouse c-MET is a heterodimer protein consisting of alpha and beta chains and contains 915 amino acid residue subunits, due to glycosylation and cleavage migrates as an approximately 45kDa (alpha chain) and 85kDa (beta chain) bands protein
Storage Instruction	Can be stored in working aliquots at 2°C-8°C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene Symbol	Met Uniprot ID= MET_MOUSE "
Immunogen Region	GLU25-ASN929

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--