

Human SIRP gamma protein (Recombinant) (STJP001879)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 M filtered PBS solution, pH7.2, 5% Trehalose.
Molecular Weight	Recombinant Human SIRP gamma (SIRP-Beta eta-2) is a protein consisting of 176 amino acid residue subunits, due to glycosylation migrates as an approximately as 45kDa band protein on reduced SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	55423
Gene Symbol	SIRPG
UniProt ID	SIRPG_HUMAN
Immunogen Region	GLU29-SER364

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--