

Mouse EGFL7 protein (Recombinant) (STJP001881)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	Insect cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 Mu m filtered PBS solution pH 7.4.
Molecular Weight	Recombinant Mouse EGFL7 is a monomer consisting of 273 amino acid residue, and migrates as an approximately 33 kDa protein under reducing conditions in SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	353156
Gene Symbol	Egfl7
UniProt ID	EGFL7_MOUSE
Immunogen Region	C-Term

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--