

Human CX3CL1/Fractalkine protein (Recombinant) (C-6His) (STJP001892)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	Pichia

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6376
Gene Symbol	CX3CL1
UniProt ID	X3CL1_HUMAN
Immunogen Region	Gln25-Gly100
Immunogen Sequence	QHHGVTKCNITCSKMTSKIP VALLIHQQNQASC GKRAII LETRQHRLFCADPK EQWVKD AMQHLD RQAAALTRNG

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

