

Human IL-34 protein (Recombinant) (C-6His) (STJP001893)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	146433
Gene Symbol	IL34
UniProt ID	IL34_HUMAN
Immunogen Region	Asn21-Pro242
Immunogen Sequence	NEPLEMWPLTQNEECTVTGF LRDKLQYRSRLQYMKHYFPI NYKISVPYEGVFRIANVTRL QRAQVSERELRYLWVLVSL S ATESVQDVLLEGHPSWKYLQ EVETLLLNQQGLTDVEVSP KVESVLSLLNAPGPNLKLVR PKALLDNCFRVMELLYCSCC KQSSVLNWQDCEVPSPQSCS PEPSLQYAATQLYPPPPWSP SSPPHSTGSRPVRAQGEGL LP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

