

## human Ephrin-A3/EFNA3 protein (Recombinant) (C-6His) (STJP001894)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">1944</a>
Gene Symbol	<a href="#">EFNA3</a>
UniProt ID	<a href="#">EFNA3_HUMAN</a>
Immunogen Region	Met 1-Ser 213
Immunogen Sequence	MAAAPLLLLLLVLPVLLPL LAQGGGALGNRHAVYWNSS NQHLRREGYTVQVNVNDYLD IYCPHYNSSGVGPGAGPGPG GGAEQYVLYMVS RNGYRTC N ASQGFKRWECNRPHAPHSP I KFSEKFQRYSAFSLGYEFHA GHEY Y Y I STPT HNLHWKCLR MKVFVCCASTSHSGEKPVPT LPQFTMGPNVKINVLEDFEG ENPQVPKLEKSI S

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

