

Human FGF-7/HBGF-7 protein (Recombinant) (C-6His) (STJP001897)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2252
Gene Symbol	FGF7
UniProt ID	FGF7_HUMAN
Immunogen Region	Cys32-Thr194
Immunogen Sequence	CNDMTPEQMATNVNCSSPER HTRSVDYMEGGDIRVRLFC RTQWYLRIKRGKVKGTQEM KNNYNIMEIRTVAVGIVAIG VGESEFYLAMNKEGKLYAKK ECNEDCNFKELILENHNTY ASAKWTHNGGEMFVALNQKG IPVRGKKTKEQKTAHFLPM AIT

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

