

## Human HMGB2 protein (Recombinant) (C-6His) (STJP001898)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">3148</a>
Gene Symbol	<a href="#">HMGB2</a>
UniProt ID	<a href="#">HMGB2_HUMAN</a>
Immunogen Region	Met1-Glu209
Immunogen Sequence	MGKGDPNKPRGKMSSYAFFV QTCREEHKKKHPDSSVNF AE FSKKCSERWKTMSAKEKSKF EDMAKSDKARYDREMKNYVP PKGDKKGGKKKDPNAPKRPPS AFFLFCSEHRPKIKSEHPGL SIGDTAKKLGEMWSEQSAKD KQPYEQKAAKLKEKEYEKDIA AYRAKGKSEAGKKGPGRPTG SKKKNEPEDEEEEEEEEEDEDED EEEEEDEE

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

