

## Active Mouse PDGF-BB protein (Recombinant) (C-6His) (STJP001923)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of 20mMNaAc-Hac pH 4.5
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">18591</a>
Gene Symbol	<a href="#">Pdgfb</a>
UniProt ID	<a href="#">PDGFB_MOUSE</a>
Immunogen Region	Ser82-Thr190
Immunogen Sequence	SLGSLAAEPAVIAECKTRT EVFQISRNLIDRTNANFLVW PPCVEVQRCSGCCNNRNVC RASQVMRPVQVRKIEIVRK KPIFKKATVLTLEDHLACKCE TIVTPRPVT

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

