

Human Chorionic somatomammotropin hormone 1/Choriomammotropin/CSH1 protein (Recombinant) (C-6His) (STJP001924)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1442
Gene Symbol	CSH1
UniProt ID	CSH1_HUMAN
Immunogen Region	Val27-Phe217
Immunogen Sequence	VQTVPLSRLFDHAMLQAHRA HQLAIDTYQEFEEYIPKDQ KYSFLHDSQTSFCFSDSIPT PSNMEETQQKSNLELLRISL LLIESWLEPVRFRLSMFANN LYYDTSDDYHLLKDLEEG IQLTMGRLEDGSRRTGQILK QTYSKFDTNSHNHDAALLKKNY GLLYCFRKMDKVVETFLRMV QCRSVEGSCGF

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

