

Human Cerebral dopamine neurotrophic factor/CDNF protein (Recombinant) (C-6His) (STJP001925)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	441549
Gene Symbol	CDNF
UniProt ID	CDNF_HUMAN
Immunogen Region	Gln27-Leu187
Immunogen Sequence	QEAGGRPGADCEVCKEFLNR FYKSLIDRGVNFSLDTIEKE LISFCLDTKGKENRLCYLLG ATKDAATKILSEVTRPMSVH MPAMKICEKLLKLD SQICEL KYEKTLDLASVDLRKMRVAE LKQILHSWGEECRACA EKD YVNL IQELAPKYAATHPKTE L

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

