

Human PDGF-BB protein (Recombinant) (No-Tag) (STJP001936)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	Pichia

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of 20mM NaAc, pH 4.5.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5155
Gene Symbol	PDGFB
UniProt ID	PDGFB_HUMAN
Immunogen Region	Ser82-Thr190
Immunogen Sequence	SLGSLTIAEPAMIAECKTRT EVFEISRRLLIDRTNANFLVW PPCVEVQRCSGCCNNRNVQC RPTQVQLRPVQVRKIEIVRK KPIFKKATVTLEDHLACKCE TVAAARPVT

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

