

Human PDGF-AA protein (Recombinant) (No-Tag) (STJP001943)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	Pichia

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mg/ml filtered solution of 20 mM NaAc, pH 4.5
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5154
Gene Symbol	PDGFA
UniProt ID	PDGFA_HUMAN
Immunogen Region	Ser87-Thr211
Immunogen Sequence	SIEEAVPAVCKTRTVIYEIP RSQVDPTSANFLIWPPCVEV KRCTGCCNTSSVKCQPSRVH HRSVKVAKVEYVRKKPKLKE VQVRLEEHLECACATSLNP DYREEDTGRPRESGKKRKRK RLKPT

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

