

Human CCL21/6Ckine protein (Recombinant) (No-Tag) (STJP001949)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	Pichia

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6366
Gene Symbol	CCL21
UniProt ID	CCL21_HUMAN
Immunogen Region	Ser24-Pro134
Immunogen Sequence	SDGGAQDCCLKYSQRKIPAK VVRSYRKQEPSLGC SIPAIL FLPRKRSQAELCADPKELWV QQLMQHLDKTPSPQKPAQGC RKDRGASKTGKKGKGSKGCK RTERSQTPKGP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

