

Human MYDGF protein (Recombinant) (C-His) (STJP001958)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 M µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	56005
Gene Symbol	MYDGF
UniProt ID	MYDGF_HUMAN
Immunogen Region	Val32-Leu173
Immunogen Sequence	VSEPTTVAFDVRPGGVVHSF SHNVGPGDKYTCMFTYASQG GTNEQWQMSLGTSEDHQHFT CTIWRPQGKSYLYFTQFKAE VRGAEIEYAMAYSKAAFERE SDVPLKTEEFVTKTAVAGR PGAFKAELSKLVIVAKASRT EL

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

