

## Human Carbonic anhydrase 5A/CA5A protein (Recombinant) (C-His) (STJP001959)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of 50mM NaAc, 50mM NaCl, 0.05% Brij-35, pH 5.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">763</a>
Gene Symbol	<a href="#">CA5A</a>
UniProt ID	<a href="#">CAH5A_HUMAN</a>
Immunogen Region	Ala40-Ser305
Immunogen Sequence	AWQTSNNTLHPLWTPVSVVPP GGTRQSPINIQWRDSVYDPQ LKPLRVSYEAASCLYIWNTG YLFQVEFDDATEASGISGGP LENHYRLKQFHFWGAVNEG GSEHTVDGHAYPAELHLVHW NSVKYQNYKEAVVGENGLAV IGVFLKLGAAHQTLQRLVDI LPEIKHKDARAAMRPFDPST LLPTCWDYWTYAGSLTTPPL TESVTWIIQKEPVEVAPSQL SAFRLLLSALGEEEEKMMV

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

