

Human Beta-NGF/NGF/NGFB protein (Recombinant) (C-His) (STJP001960)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4803
Gene Symbol	NGF
UniProt ID	NGF_HUMAN
Immunogen Region	Glu19-Ala241
Immunogen Sequence	EPHSESNVPAGHTIPQAHWT KLQHSLDTALRRARSAPAAA IAARVAGQTRNITVDPRLFK KRRLRSPRVLFSTQPPREAA DTQDLDFEVGGAAPFNRTHR SKRSSHPIFHRGEFSVCDS VSVWVGDKTTATDIKQKEVM VLGEVNINNSVFKQYFFETK CRDPNPVDSGCRGIDSKHWN SYCTTTHTFVKALTMGKQA AWRFIRIDTACVCVLSRKAV RRA

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

