

Human CER1/DAND4 protein (Recombinant) (C-His) (STJP001961)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of 25mM NaAc pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9350
Gene Symbol	CER1
UniProt ID	CER1_HUMAN
Immunogen Region	Thr18-Ala267
Immunogen Sequence	TRHQDGRQNQSSLSPVLLPR NQRELPTGNHEEAEEKPDLF VAVPHLVATSPAGEGQRQRE KMLSRFGRFWKKPEREMHPS RDSDSEPFPPGTQSLIQPID GMKMEKSPLREEAKKFWHHF MFRKTPASQGVILPIKSHEV HWETCRTVPFSQTITHEGCE KVVVQNNLCFGKCGSVHFPG AAQHSHTSCSHCLPAKFTTM HLPLNCTELSSVIKVVMLVE ECQCKVKTEHEDGHILHAG

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

